

Real time monitoring for gas phase H/D exchange reaction in H₂O and D₂O by Using a Compact High-Resolution Multi-Turn Mass Spectrometer

AN-MS-002

Introduction

The InfiTOF™ is a compact time-of-flight mass spectrometer which has unique multi-turn ion optics.

The InfiTOF has features below.

1. Able to detect small elements such as Hydrogen atom (H).
2. Allow resolving power selection based on the chosen path length for the ions.
3. High resolving power achieved to separate close proximity ions such as carbon monoxide (CO) and nitrogen (N₂).
4. Able to measure long term with high m/z stability (Fluctuation < 100ppm/hour), high sensitivity and high mass accuracy.
5. Easily transportable.

The InfiTOF is suitable for real time gas reaction monitoring. In this report, we showed real time monitoring for H/D exchange reaction of light water (H₂O) and heavy water (D₂O) which is one of the simplest chemical reactions.

Experimental

Mass spectra were acquired by using the InfiTOF compact multi-turn high-resolution mass spectrometer. High-resolution mass spectra were acquired in multi-segment mode (150 turns) to give a resolving power of approximately 30,000 (FWHM) at m/z 18. Additionally, the ionizing electron voltage was set to 15 eV. Schematic drawing of experimental equipment shown in Fig.1. Between the gas chromatograph and the mass spectrometer, a vial of volume 22 mL was set up, and the two intervals as GC – vial and vial – MS were connected through deactivated fused-silica capillary tubes (3m x 0.25mm and 30cm x 0.1mm respectively). GC Helium gas flow rate was set up 1mL/min.

H₂O and D₂O were injected 1μL each in GC inlet at a one minute time lag. Introduced gases were through in the vial, and the gases were mixed in the vial were introduced into the mass spectrometer.




Figure 1. Schematic drawing of experimental equipment.




Table 1 Measurement conditions

Sample	H ₂ O (W5-4, Fisher scientific, purified by EASYpure RF, Barnstead™)
	D ₂ O (DLM-6-1000, Cambridge Isotope Laboratories, Inc.)
Mixing vial volume	22mL
Vial temperature	43°C
Injection volume	1μL each
Injected timing	H ₂ O ; Injected 30 second after started to measure
	D ₂ O ; Injected 90 second after started to measure
MS	JMS-MT3010HRGA InfiTOF™
Resolving power	Approx. 30000 (m/z 18)
Ionization Voltage	15 eV
Recording interval	1 sec/spectrum
Detector Voltage	2700V
Measurement mode	Extended mode
m/z range	m/z 17 – 22
GC	Agilent 7890B
Column (GC to Vial side)	Deactivated fused silica capillary tube (3m x 0.25mm)
Column (Vial to MS side)	Deactivated fused silica capillary tube (30cm x 0.1mm)
He flow rate	1 mL/min

Results

The extracted ion chromatogram in Fig.2 shows that each ions behavior is different. These ions were detected for 1 hour after injecting light water. H/D exchange reaction was occurred instantaneously when injecting heavy water seems to be detected HDO⁺ and D₂O⁺.

The mass spectra after light water injection and after light water and heavy water injection respectively are shown in Fig. 3. The mass spectrum that injected only light water shows H₂O⁺ and H₃O⁺ were detected mainly. The mass spectrum after both water injections showed the ions which mass difference is only 0.0015 u such as HDO⁺ /H₃O⁺ and D₂O⁺ /H₂DO⁺ were detected respectively.




Figure 2. Extracted ion chromatogram.





Figure 3. Mass spectra (Left; After injected only H₂O, Right; After injected H₂O and D₂O).

Conclusion

InfITOF can separate and identify to close proximity ions such as HDO/H₃O and D₂O/H₂DO by high mass resolving power. Therefore the InfITOF is suitable system for absolute H₂O/D₂O ratio analysis.