

A ^1H Solid State NMR Application for Pharmaceutical Development by using Ultra Fast MAS : Drug-Polymer Intermolecular Interaction on Solid Dispersions

Product used : Nuclear Magnetic Resonance(NMR), 1mm HX probe

The low oral bioavailability of a drug due to its poor aqueous solubility is a major challenge for pharmaceutical development. Solid dispersion (SD), where the amorphous drug is dispersed into the polymer matrix, is one of the useful approaches to improve the aqueous solubility. However, thermodynamically unstable nature of an amorphous drug increases its susceptibility to recrystallize upon storage, which, in turn, reduces its solubility and dissolution. Therefore, design of thermodynamically stable SD is required. One of the important points in the design is intermolecular interaction, such as hydrogen bond, between active pharmaceutical ingredient (API) and amorphous polymer. Solid state NMR, particularly on ^1H solid state NMR using ultra fast magic spinning, is a powerful tool to probe intermolecular interactions because it is able to observe hydrogen bond directly through hydrogen atoms. We report here an example of structural analysis for a carbamazepine (CBZ) / hypromellose (HPMC) SDs. ⁽¹⁾

CBZ and HPMC derivatives, including HPMC, hypromellose acetate (HPMC-A), and hypromellose succinate (HPMC-S), were spray-dried to prepare CBZ/polymer spray-dried samples (SPDs). CBZ/HPMC and CBZ/HPMC-A SPDs recrystallized within 10 days, whereas CBZ/HPMC-S SPD maintained its amorphous state for a longer period. This indicates that amorphous state of CBZ/HPMC-S SPD is more stable than the other SPDs.

Fig.2 shows ^1H single pulse MAS NMR spectra of CBZ/HPMC, CBZ/HPMC-A, CBZ/HPMC-S (CBZ loading of 33 wt%), and HPMC-S SPD. It is noteworthy that only CBZ/HPMC-S SPD exhibited an ^1H peak at a magnetic field lower than 10 ppm. The fact HPMC-S SPD did not show an ^1H peak at this field indicates that the COOH group on HPMC-S forms hydrogen bond with the C=O group of CBZ in CBZ/HPMC-S SPD.

To further investigate the state of NH_2 in CBZ, ^1H spectra through ^1H - ^{15}N Cross polarization (CP) based HSQC filter were obtained by using ^{15}N labeled CBZ. By setting a short enough second contact time $50\mu\text{s}$ (^{15}N - \rightarrow ^1H), NH_2 proton covalently bonded to ^{15}N nucleus can be observed selectively. The HSQC filter ^1H NMR spectra of CBZ/HPMC, CBZ/HPMC-A, and CBZ/HPMC-S SPDs with ^{15}N labeled CBZ loading of 33 wt% are shown in Fig. 3. A single ^1H peak was observed at 5 – 6 ppm in the spectrum of CBZ/HPMC SPD. In contrast, CBZ/HPMC-A and CBZ/HPMC-S SPDs showed a small peak shoulder at 6 – 8 ppm. This implies the formation of hydrogen bonding between the NH_2 groups of CBZ and the polymers. Furthermore, the ^1H peak shoulder was stronger in CBZ/HPMC-S than that in CBZ/HPMC-A, indicating a larger proportion of CBZ molecules in CBZ/HPMC-S formed CBZ – polymer interactions as compared to those in CBZ/HPMC-A. Interactions between CBZ and HPMC polymers are schematically summarized in Fig. 4. Solid State ^1H NMR is a potent technique to investigate drug-polymer intermolecular interactions on SDs, and this leads to understand relationship between their structures and stability at molecular level.




Fig.1 Chemical structure of (a) CBZ and (b) HPMC




Fig.2 ^1H single pulse MAS NMR spectra of CBZ/HPMC, CBZ/HPMC-A, CBZ/HPMC-S, and HPMC-S SPDs. (Reprinted with permission from ref 1. Copyright at 2018 American Chemical Society)




Fig.3 ^1H - ^{15}N CP-based HSQC filter ^1H MAS NMR spectra of CBZ/HPMC, CBZ/HPMC-A, and CBZ/HPMC-S SPDs with ^{15}N -labeled CBZ loading of 33 wt%. The solid arrow in this figure represents the newly developed peak shoulder.




Fig.4 A schematic illustration of CBZ-HPMC interactions in CBZ/HPMC SPDs. (Reprinted with permission from ref 1. Copyright at 2019 American Chemical Society)

Reference:

1) Y. Ishizuka, K. Ueda, H. Okada, J. Takeda, M. Karashima, K. Yazawa, K. Higashi, K. Kawakami, Y. Ikeda, K. Moribe, *Mol. Pharmaceutics.*, **2019**, *16*, 6, 2785–2794

Copyright © 2019 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

