

Solutions for Innovation



Gas analysis solutions

Gasanalysis

Introduction

Gas analysis is important for numerous applications including, characterization of impurities in highpurity gases, semiconductor manufacturing, gases released during battery charge and discharge cycles or from biomass, combustion, pyrolysis of polymer materials, and nuclear reactors. JEOL mass spectrometers offer gas analysis solutions to match the needs of various processes and research fields.



A comparison of JEOL gas analysis systems



GC-QMS JMS-Q1500GC



High Performance GC-TOFMS JMS-T200GC AccuTOF[™] GCx-plus

Multi-Turn TOFMS JMS-MT3010HRGA INFITOF

Mass resolution	2,000	10,000	30,000 *1	Applications
Separation of nitrogen (N ₂) and carbon monoxide (CO ₂) (m/z 28)	X	$\sqrt{\sqrt{2}}$	<i>√√√</i>	Quantitation of carbon monoxide (CO) as an impurity in a high purity gas. Tracking a reaction that generate or consume CO. (Quantitation of CO without being affected by nitrogen (N_2) background.)
Separation of nitrous oxide (N ₂ O) and carbon dioxide (CO ₂) (m/z 44)	X	X	$\sqrt{\sqrt{\sqrt{1}}}$	Simultaneous analysis of nitrous oxide (N ₂ O) and carbon dioxide (CO ₂).
Detection of hydrogen related ions				
Hydrogen molecular ion $(H_2^+; m/z 2)$	$\sqrt{\sqrt{2}}$	X	$\sqrt{\sqrt{2}}$	Quantitaion of hydrogen (H ₂) as an impurity in a high purity gas. Thermal desorption of H2 from metal. (EGA-MS, TG/DTA-MS)
Hydron (hydrogen atomic ion) (H $^+$, <i>m</i> / <i>z</i> 1)	X	X	$\sqrt{\sqrt{\sqrt{1}}}$	Tracking a reaction involving hydrogen radicals.
Transportability	\checkmark	X	$\sqrt{}$	On-site analysis in the vicinity of the gas source.
High flow rate introduction	$\sqrt{\sqrt{2}}$	\checkmark	$\sqrt{}^{*2}$	Tracking a rapid reaction. Monitoring combustion.
Prolonged continuous monitoring	$\sqrt{\sqrt{\sqrt{1}}}$	\checkmark	$\sqrt{\sqrt{\sqrt{1}}}$	QA/QC of a high purity gas. Prolonged monitoring of a reaction process.
EGA (evolved gas analysis) (with multi-functional pyrolyzer)	$\sqrt{\sqrt{\sqrt{1}}}$	$\sqrt{\sqrt{\sqrt{1}}}$	* "3	Testing thermophysical properties of an organic or inorganic material.
TG/DTA-MS	$\sqrt{\sqrt{\sqrt{1}}}$	$\sqrt{\sqrt{\sqrt{2}}}$	$\sqrt{\sqrt{\sqrt{1}}}$	
GC/MS with a packed column	$\sqrt{\sqrt{\sqrt{1}}}$	X	X	Analysis of impurities in a high-purity gas.
GC/MS with a capillary column	$\sqrt{\sqrt{2}}$	$\sqrt{\sqrt{\sqrt{1}}}$	* "3	

Legend

√√√: Excellent

√√: Good

√: Possible

X: Not possible

*1 Mass resolution is variable according to the flight cycles.
*2 Optional high volume rotary pump or scroll pump is required.

*3 Please inquire.

Quadrupole Mass Spectrometer

JMS-Q1500GC is the fifth generation of the JEOL's GC-QMS system. It employs large hyperbolic quadrupole and high capacity vacuum pumping system, enabling stable quantitative analysis of hydrogen. Since it achieves high sensitivity under high carrier gas flow rate, it is ideal for GC/MS analysis with packed column, on-line connection with thermogravimetry (TG) and temperature programmed desorption (TPD).



The quadrupole GC-MS system can analyze gas, liquid and solid samples and detect volatile compounds in the specimen and/ or generated from the specimen.

A variety of applications are available to analyze inorganic gases, organic compounds in solution and additives included in materials such as resin.

Analysis example Prolonged

Prolonged monitoring of a standard gas

Sample : Standard gas; 10 ppm each of H_2 , O_2 , N_2 , CH_4 , Ar, and CO_2 in He

m/z 28 : CO⁺⁺/N₂⁺⁺ \Rightarrow Both species contribute to the peak height because they are not separated at the mass resolution of a quadrupole mass spectrometer.

m/z 18 : H_2O^{+-} \Rightarrow Water in the vacuum chamber and tubing.



Maximum time for one measurement is 24 hours. Monitoring longer than 24 hours is possible by automatically repeating the measurement.



Analysis of trace hydrogen in a gas mixture by using packed-column gas chromatography

Chromatography maximizes the utility of the system for the analysis of gas mixtures. JMS-Q1500GC with high capacity vacuum pumping system can easily handle high carrier gas flow late of packed GC columns, which are essential for the analysis of trace components in gas mixture.

Analysis condition

MS Condition	
Ionization mode	EI (17 eV, 30 μΑ)
Mass range	<i>m/z</i> 2-200
Ion source temperature	200 °C
GC interface temperature	200 °C
GC condition	
Injection mode	Direct to packed column
Inlet temperature	200 °C
Column	A proprietary packed column
Oven temperature	40 °C (5 min) \rightarrow 8 °C /min \rightarrow 220 °C
Sample injection	6-port valve Loop size:250 μL

TIC Chromatogram (scanning mode)



Extracted Ion Chromatogram (scanning mode) of H_2^{+} (m/z 2)









Scan measurement for identifying unknown compounds and selected ion monitoring (SIM) for target compound quantitation can be performed simultaneously. (SCAN/SIM function)

High Performance Gas Chromatograph Time-of-Flight Mass Spectrometer

JMS-T200GC AccuTOFTM GCx-plus is the 5th generation of JEOL's successful AccuTOFTMGC series of GC-HRTOFMS systems. With the patented ion transfer system that achieves high ion transmission without using radio-frequency-driven ion guide, high sensitivity is achieved even for low molecular weight gas components. It is challenging to achieve high mass-resolving power at low *m*/*z* region on a short flight path time-of-flight mass analyzer. JMS-T200GC AccuTOFTM GCx-plus employs high speed (4 gigasamples/s) digitizer and newly developed high speed ion detector, achieving high mass-resolution and excellent mass spectral peak shape even in low *m*/*z* region.

Several ion sources are available including Electron Ionization (EI, standard), and optional Chemical Ionization (CI), Field Ionization (FI), and Photoionization (PI) sources. This offers the widest variety of analyses for volatile compounds.



The High Performance Gas Chromatograph Time-of-Flight Mass Spectrometer can acquire mass spectra with high sensitivity and high speed while providing high mass resolving power and mass accuracy for samples introduced by gas chromatography.









Separation of isobaric components at m/z 28



Even at m/z as low as 28, mass resolution sufficient to resolve CO and N₂, and excellent mass spectral peak shapes are achieved.

Analysis example

Evolved gas analysis (EGA) of nylon 66 (PA66) by using a multi-functional pyrolyzer

When synthetic polymers are heated to high temperature, low molecular weight gas evolve. With JMS-T200GC AccuTOF[™] GCx-plus, elemental composition of each evolved gas component can be elucidated through accurate mass measurement. High mass-resolution extracted ion chromatograms facilitate to monitor time course change of each evolved gas component.

Analysis conditions

Instrument	JM S-T200GC "AccuT OF™ GCx"
Py Conditions	
Pyrolysis Furnace Temperature	60°C ⇒ 20°C/m in ⇒ 700°C (a m inute)
Py-GC-ITF Temperature	350 °C
GC Conditions	
Inlet temperature	350°C
Injection mode	Split 50 :1
Columns	Deactivated fused silica tube (5 m x 0.25 m m)
Oven Temperature Program	350°C (33 m in)
Carrier Gas flow	1 m L/m in (He, Constant flow mode)
MS Conditions	
Ionization method	EI(+); 25 V, 200 μA
Interface Temperature	350°C
Ion source Temperature	250°C
Spectrum recording intervals	2.5 Hz (0.4 sec / Spectrum)
<i>m/z</i> range	10~800
Drift Compensation	m/z 28.0056 (N2 ⁺)
Sample Weight	0.28 mg
	L H Ojn
	PA66 (Nylon 66)
	Measured Value Elemental Composition Calculated Mass [mDa]

Mass spectrum



18.0104

18.0344

27.9950

H₂O⁺⁻

NH₄+

C0+-

CH₂N

18.0100 0.4

0.6

18.0338 0.6

27.9944

EGA total ion current thermogram and mass spectrum



High mass-resolution extracted ion chromatograms



Multi-Turn Time-of-Flight Mass Spectrometer

Multi-Turn Time-of-Flight Technology is a breakthrough technology for compact high-resolution mass spectrometers. This compact TOF-MS is designed in response to cutting edge research needs for in-situ realtime analysis of small gas samples. The system offers high mass resolving power for the analysis of advanced materials. With a wide mass range (m/z 1-1000), the system is ideal for the research related to next-generation energy sources and hydrogen fuel.





Detection of hydrogen

With conventional high-resolution TOFMS, detecting hydrogen ions is highly challenging. In contrast, the INFITOF can detect hydrogen ions (H⁺ at m/z 1 and H₂⁺⁺ at m/z 2) without sacrificing the ability to detect larger species (up to m/z 1000).



Analysis example

Continuous monitoring of low concentration gas mixture containing hydrogen; Resolution \approx 5,000 (4 hours)

Stable signal over a 4-hour period for gases including hydrogen Measurement Condition: Ionization: EI, Ionization voltage: 25 eV, Analyzer : 14 turns (cycles), $R \ge 5,000 \ (m/z \ 28)$







Analysis example

Continuous monitoring of high concentration gas containing hydrogen and N₂O; Resolution \approx 12,000

Stable analysis of high concentration gas in the % order. Complete separation of N_2O (*m/z* 44.00106) and CO_2 (*m/z* 43.98983) is highly effective for the analysis of combustion products.

Measurement condition : Ionization : EI, Ionization voltage : 25 eV, Analyzer : 30 turns (cycles), R ≥ 12,000 (m/z 28)



TG/DTA-MS System

TG/DTA (Thermogravimetry / Differential Thermal Analysis) system records the mass changes of and endothermic/exothermic heats to/from the specimen while the temperature of the specimen is increased. By combining the TG/DTA system with a mass spectrometer, components evolved from the specimen can be detected and identified.



Pyrolyzer-GC-MS System

A Pyrolyzer-GC-MS System, which combines a Pyrolyzer with a GC-MS, is ideal for the structural analysis of synthetic polymers by pyrolisis and thermal extraction of semi-volatile components in materials. Pyrolysis GC/MS is widely used for materials analysis.









Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

JEOL DEOL Ltd.

3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577 www.jeol.com ISO 9001 • ISO 14001 Certified

AUSTRALIA & NEW ZEALAND /JEOL (AUSTRALASIA) Pty/Ltd. Suite 1, L2 18 Aquatic Drive - Frenchs Forest NSW 2086 Australia • BELGIUM /JEOL (EUROPE) B.V. Planet II, Gebouw B Leuvensesteenweg 542, B-1930 Zaventem Belgium /BOL Challo Brazil Instrumentos Clentificos Ltda. Av. Jabaquara, 2965 6° andar conjunto 52; 0404-500 Sao Paulo, SP Brazil - CANADA /JEOL (JERVICE BUREAU 3rd FI. Nile Center Bidg, Navai Street, Doki, (Jerone) FRANCE /JEOL (JERVICE BUREAU 3rd FI. Nile Center Bidg, Navai Street, Doki, (Jerone) - Cegnton / Allee de Giverny 77290. Croixey-uc-Seine, France • GERMANY /GmUH GULA Anger 20 85365 Friesing. Germany • GREAT BRITIN & IRELAND /JEOL (JUK.) J. JEOL House. Silver Canada / USA /JEOL // HALYSIA /JEOL (MALYSIA) SDI. A. Palazzo Pacinotti - Milano 3 City, Via Ludovico il Moro, 6/A 20080 Basiglio(MI) Italy • KOREA /JEOL KOREA LTD. Dongwoo Bldg. 7F, 1443, Yangjae Daero, Gangdong-Gu, Segadon Benic Juarez, C.P. (38) Block /A Level 5, Kelana Business Center, 97, Jalan S5 7/2, Kelana Jaya, 47301 Petaling Jaya, Selangor, Malaysia /JEOL (MALYSIA) SDI. MEXICO JACO • CATAR /Mannai Trading Company / ULL. ALL Ematid Complex, Salva Road P.O.Box 76, Dong, Qatar • NBLAYSIA /JEOL (HUSI) LEC Krasnoprotelarskaya Street, 16, Bld. 2, 127473, Moscow, Russian Federation • SCANDINAVIA /SWEDEN JEOL (Nordic) AB Hammarbacken 6A, Box 716, 191 27 Sollentuna Sweden • SINGAPORE /JEOL ASIA PTE.LTD. 2 Corporation Place Singapore 61494 • TAIWA /JIE DONG CO., LTD. 7F, 112, Chung Hsiao East Road, Section 1, Taipel, Taiwan 10023 (R.O.C.) • THE NETHERLANDS /JEOL (EUROPE) B.V. Lireweg