

Introduction

Context

In the world of microscopic studies, the **identification of the spatial distribution** of certain mineral phases and inclusions **remains complex**. Indeed, variations in **chemical composition** on a scale that can sometimes be very small remain **difficult to observe** with a simple optical microscope and even with the use of a scanning electron microscope (SEM). The **identification** and the **precise localization** of these differences of chemical nature remain however **essential for the good understanding** and the interpretation of the sample, and those, even on a **very small scale** (going to the order of the hundred or even of the ten nanometer).

Case of study

In this study, a geological thin section of a **meta-arkose in a serpentinite** was chosen. This sample comes from the Bureau de Recherches Géologiques et Minières (BRGM). This type of metamorphic rock comes from the interaction between sea water and peridotites, which are the rocks that make up the earth's mantle. Within the meta-arkose rock is a **pyrite surrounded by different mineral phases** that appear to have **variable chemical composition**. This area was therefore selected for the **identification and distribution of the different elemental assemblages**.


Methodology

In order to better answer the questions of chemical identification and spatial distribution, a **set of multi-approach and complementary studies** was conducted between JEOL and HORIBA companies. It includes the use of a **SEM IT800HL** of the JEOL brand to which were added three other tools, namely an **Energy Dispersive X-Ray Spectroscopy (EDS)** and a **Soft X-Ray Emission Spectrometer (SXES)** of the JEOL brand as well as a versatile hyperspectral **cathodoluminescence (CL) F-CLUE** of the HORIBA brand. In addition, the sample selected for this study was analyzed by **RAMAN** with the LabRAM Soleil™ Raman Microscope from HORIBA.

Main results

SEM - BSE image



Backscattered electrons (BSE) are high-energy electrons that are produced by the elastic scattering of the primary beam electrons of the SEM with the atom nuclei. The yield of BSE, that is the ratio of the number of emitted BSE and the amount of primary beam electrons, **depends on the atomic number**: the higher the atomic number, or the heavier the element, the brighter the contrast.



In this study, BSE imaging allows to see the **difference of chemical composition** on the observed area. Thus, this study area showing **important contrasts**, sign of **important chemical variation**, is selecting for additional on complementary analysis.

Raman mapping


Raman photons are emitted when the sample is illuminated by a laser source (UV-visible-IR) through an inelastic light scattering phenomenon. The **gain or loss of energy** of the inelastic photons emitted compared to the incident photons is translated on the **Raman spectra by a frequency shift**. The Raman bands observed at a frequency shift correspond to the **energy gap between the vibrational levels of the analyzed molecule**.



In this study, the RAMAN allows the **identification of the main minerals** composing and surrounding the pyrite as well as their **general spatial distribution**.

Nano GPS NavYX

This GPS allows to spot the same points of interest with a **high degree of accuracy** and **simplicity** using different microscopes. It is based on **patterned tags** fixed to the samples. An image of a part of the tag can be automatically **converted to absolute coordinates** and **angular orientation** which will provide a **correspondence** between the **sample** and the coordinates of the **moving stage**.



With NavYX, maps from different microscopic techniques are **co-localized** and overlaid to create **composite maps**. Thus, we can perform **correlative measurements** on the same observation area.

